Angio cpet - with automatic stand



Modern ergometer with multifunctional applications







Highlights

Easy to operate

For Lode products this means:

- easy to connect
- easy to move around
- easy user interface

Reliable and reproducible stress tests

The experience of professionals who calibrate many ergometers shows that the Lode ergometers are the most reliable across the complete workload and rpm range and still within specifications even after many years of intensive use.

High standards

Lode is a socially and environmentally responsible company. All Lode products are RoHS/WEE compliant and Lode is ISO 9001:2003, ISO 13485:2008 and FDA 510K certified. All medical products comply to MDD 93/42/EEC, incl. IEC 60601-1.

Additional features with PCU

Besides the possibility to program 24 protocols easily, this control unit offers the following features:

- better monitoring because of the additional and larger display
- a perfect combination with BPM
- possibility to measure SpO2

Multifunctional

The ergometer can be used for both arm and leg ergometry.





Angio cpet - with automatic stand



Modern ergometer with multifunctional applications

With the stand for the Angio, the ergometer can move up and down over a range of 140cm. With the stand it is easier to move your arm Ergometry setting to another room, because it can stand free. For heavy use during arm-ergometry it is possible to fixate the stand to the wall with the supplied parts. Together with the automatic stand the Angio becomes an ergometric unit that can be used for both arm and leg Ergometry (this setting is standard supplied with handcurbles). Its compact design makes it universally applicable for ergometric studies in those sectors in which standard ergometry cannot be used. The Angio operates independent of pedaling speed in the range of 7 - 1000 watt. The Angio cpet is standard supplied with a communication module and can therefor be easily controlled by all known stress ECG and pulmonary devices in the world. The workload, rpm and time can be readout from the 3,5" colour display.

Features



Small adjustment steps

The workload of the Lode ergometers is adjustable in steps of only 1 watt. Depending **//211** on your wishes, the test operator or the test subject can adjust the workload. The steps of 1 watt are possible in the manual mode as well as within protocols.



Extreme low start-up load

The extreme low start-up load of 7 watts and the adjustability in small steps of 1 watt make **1211** this ergometer perfectly suitable for many different applications. The standard control unit shows multiple ergometry parameters and you can determine your specific default setting and start-up menu.



Accurate over a long period of time

The Lode ergometers are supplied with an electro-magnetic braking mechanism of Lanooy (eddy current). The biggest advantage of this braking system compared to a friction braking system is the absolute accuracy and the accuracy over time. Moreover, friction braking systems have more wearing parts.



LEM compatible

This product can be used with Lode Ergometry Manager (LEM) software to manage data and to apply specific protocols when a Communication card or Network card is present



LCRM compatible

This product can be used with Lode Cardiac Rehabilitation Manager software (LCRM)



Customer specific display setting

Display settings are adjustable according to your specific requirements: each individual has its specific wishes about the parameters to be displayed. This can easily be adjusted with the Lode ergometers.



Versatile Interfacing

Various interface protocols guarantee perfect communication with all commonly known stress ECG and spirometry equipment



Angio cpet - with automatic stand



Modern ergometer with multifunctional applications

Angio cpet - with automatic stand can a.o be extended with the following options:

Control Unit with touch screen 7" for ergometer Multifunctionality



Partnumber: 945834

Programmable Control Unit with 7" Touchscreen Programmable



Partnumber: 945835

Heart rate large distance

Heart rate in beats per minute



Partnumber: 945830

SpO2 for control unit with touch panel (bicycle)

Saturation and heart rate



Partnumber: 945823

Stabilisation leg for Angio with Automatic Stand

Optimal safety and stability



Partnumber: 967815

Extra wide access for wheelchairs

Easy access



Partnumber: 967805

Pedal shoes (pair)

Extra stability during cycling



Partnumber: 917803

Pedal shoes pediatric (pair)

Pedal shoes for childen



Partnumber: 917833

diatric Pedal shoes extra large (pair)

For large feet sizes



Partnumber: 917834

Adjustable cranks

Optimal force application



Partnumber: 928804

Handgrip fixation set

Handgrip fixation set



Partnumber: 917822

Handgrip fixation set pediatric

Versatile ergometry



Partnumber: 917823

Electric adjustable chair for arm ergometer

Comfortable seating position in front of the



Partnumber: 917813

Network Module

Easy networking with LEM and LCRM



Partnumber: 945851

RS232 cable

Easy connection



Partnumber: 930911



Angio cpet - with automatic stand



Modern ergometer with multifunctional applications

Specifications

Workload		
Minimum load	7 W	
Maximum peak load	1000 W	
Minimum load increments	1 W	
Maximum continuous load	750 W	
Hyperbolic workload control	~	
Maximum rpm independent constant load	150 rpm	
Minimum rpm independent constant load	30 rpm	
Electromagnetic "eddy current" braking system	~	
Dynamic calibration	~	
Accuracy		
Workload accuracy from 7 to 100 W	3 W	
Workload accuracy from 100 to 500 W	3 %	
Workload accuracy from 500 to 1000 W	5 %	
Comfort		
Adjustability range height	1400 mm	
User Interface		
Readout Distance	~	
Readout RPM	~	
Readout Heartrate	~	
Readout target HR	~	
Readout Energy	~	
Readout Torque	~	
Read out Time	~	
Readout Power	~	
Set Display	~	
Set Resistance	~	
Set P-Slope	~	
Set Mode	~	
Manual operation mode	~	
Preset protocol operation mode	~	
Terminal operation mode	~	
External control unit	~	
Selfdesigned protocol operation mode	~	
Dimensions		
Product length (cm)	84 cm	33.1 inch
Product width (cm)	82 cm	32.3 inch
Product height	223 cm	87.8 inch

Power requirements

VAC	100 - 240 V
Phases	1
Frequency	50/60 Hz
Power consumption	250 W
Power cord IEC 60320 C13 with CEE 7/7 plug	~
Power cord NEMA	×
Standards & Safety	
IEC 60601-1:2012	~
ISO 13485:2003 compliant	~
ISO 9001:2008 compliant	~
Certification	
CE class Im according to MDD93/42/EEC - pending	~
CTüVus according to NRTL - pending	~
CB according to IECEE CB - pending	~

Order info

Product weight

Width between stand

Partnumber: 967904

*Specifications are subject to change without notice.



Distributed by

Costa Rica

Costa Rica Tel: Lode B.V. Zernikepark 16 9747 AN Groningen The Netherlands Tel: +31 50 5712811 Fax: +31 50 5716746 E-mail: ask@lode.nl

Internet: www.lode.nl

209.4 lbs

95 kg

67 cm